



PERFORMANCE QUALITY SAFETY



E008: Auto-disable syringe for fixed dose immunization

PQS code: E008/051
Description: AD syringe 0.5ml
Manufacturer's reference: JMI AD syringe 0.5ml
Manufactured in: Bangladesh
Company: JMI Syringes & Medical Devices Ltd.
Address: Unique Heights, Level-7
117, Kazi Nazrul Islam Avenue
Ramna, Dhaka-1217
Bangladesh
Telephone: +880-8170681 to 5
Email: info@jmigroup-bd.com
Web address: <http://www.jmisyringe.com/>

Specifications

Graduations:	0.5ml	Quality standard:	ISO 13485,
Syringe material(s):	Polypropylene	Markings:	CE mark
Fixed needle size:	23G x 1" (0.60 x 25 mm)	Pieces per carton:	2400
AD mechanism:	Plunger lock	Volume per carton (m3):	0.133 m3
AD location:	Upon completion of injection	Weight per carton (kg):	15.6 kg
Number of components:	3 pieces	Minimum order:	192,000
Other needle options:	No	Incoterms:	FOB
Primary packaging:	Blister pack	Prices per unit:	Upon request to manufacturer
Year base price:	2013		

Comments:

Current PQS status: pre-qualified ; 29 May 2013

Valid until : May 2026

Note: If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased